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TO AVOID DANGER OF SUFFOCATION, KEEP THIS PLASTIC BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE THIS BAG IN CRIBS, BEDS, CARRIAGES OR PLAYPENS. THIS BAG IS NOT A TOY. Nebezpečí udušení.

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Vyrobeno: Guangdong Transtek Medical Electronics Co., Ltd.
Zone A, No.105 ,Dongli Road, Torch Development District
Zhongshan,528437,Guangdong,China

Notifikovaná osoba pro Evropu

MDSS - Medical Device Safety Service GmbH; Schiffgraben 41, 30175 Hannover, Germany

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eta 4297

INSTRUCTIONS FOR USE

Dear customer, thank you for purchasing our product. Please read the operating instructions carefully before putting the appliance into operation and keep these instructions including the warranty, the receipt and, if possible, the box with the internal packing.

I. SAFETY WARNING

- Consider the instructions for use as a part of the appliance and pass them on to any other user of the appliance.
- **The appliance is intended for use in households only and for similar use!**
It is not intended for use in health-care facilities and for commercial use! The values measured are intended for your information and monitoring the trends only and must not in any way replace a visit to the physician and having the blood pressure checked under qualified supervision. Proper blood pressure measurement must be performed by a qualified person (physician, nurse etc.) and in the context of an overall health condition of the person.
- This appliance must not be used by children. This appliance can be used by persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision. Children younger than 8 years must be kept out of reach of the appliance and its power cord.
- To power the appliance, only use the adapter included.
- If the adaptor of the appliance is corrupted, it has to be replaced by the manufacturer, its service technician or a similarly qualified person so as to prevent dangerous situations.
- Never use the appliance if its adaptor or plug is damaged, if it does not work properly, if it fell down and was damaged. In this case take the appliance to a special service to check its safety and proper function.
- Before replacement of accessories or accessible parts, before assembly and disassembly, before cleaning or maintenance, switch off the appliance and unplug adaptor from the mains!
- Always unplug the appliance from charger if you leave it unattended.
- The appliance with the plugged adapter must not be used at places where it could fall into a bath tub, a wash basin or a swimming pool. If the appliance fell to water anyway, do not take it out! First of all, unplug the adapter from the socket and only then remove the appliance. In this case take the appliance to a special service to check its safety and proper function.

- To power the appliance, only apply safe small voltage in accordance with the type plate.
- Do not put or take out the plug of the adapter to or from a socket with wet hands!
- The adaptor supplied with this appliance is intended for this appliance only, do not use it for other purposes. Also, use only the supplied adaptor with this appliance, do not use any other adaptor for charging (e.g. from another appliance).
- Never immerse the device or adapter in water or other liquids (nor partially), do not expose them to moist and do not wash them under the tap!
- The power cord of the adapter must not be damaged with sharp or hot objects, open fire and it must not be sunk into water.
- Please note that the data obtained from the device are only approximate and inaccurate in comparison with the results of real medical analyses. Precise assessment can be made only by a physician using appropriate instruments.
- Do not use rechargeable accumulators.
- If the battery is leaking, remove it immediately; it can damage the device.
- Dispose of a discharged batteries in a suitable way (see par. **VII. ENVIRONMENTAL PROTECTION**).
- Keep batteries and appliance out of reach of children and incapacitated people. A person, who would swallow the batteries, must seek medical help immediately.
- All text in other languages and pictures shown on the packaging, or the product itself, are translated and explained at the end of this language mutation.
- **WARNING:** There is a risk of injury in case of incorrect use (not in accordance with the manual) of the appliance.
- The producer is not responsible for damage caused by improper use of the appliance and its guarantee does not cover the appliance in the case of failure to comply with the safety warnings.
- This device is intended for adult use in homes only.
- The device is not suitable for use on neonatal patients, pregnant women, patients with implanted, electronic devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral, arterial disease and patients undergoing intravascular therapy or arterio-venous shunt or people who received a mastectomy. Please consult your doctor prior to using the unit if you suffer from illnesses.
- The device is not suitable for measuring the blood pressure of children. Ask your doctor before using it on older children.
- The device is not intended for patient transport outside a healthcare facility.
- The device is not intended for public use.
- This device is intended for no-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Do not begin or end medical treatment without asking a physician for treatment advice.
- If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your physician.
- Do not take any therapeutic measures on the basis of a self measurement. Never alter the dose of a medicine prescribed by a doctor. Consult your doctor if you have any question about your blood pressure.
- Don't kink the connection tube during use, otherwise, the cuff pressure may continuously increase which can prevent blood flow and result in harmful injury to the PATIENT.

- When using this device, please pay attention to the following situation which may interrupt blood flow and influence blood circulation of the patient, thus cause harmful injury to the patient: connection tubing kinking too frequent and consecutive multiple measurements; the application of the cuff and its pressurization on any arm where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present; inflating the cuff on the side of a mastectomy.
- **Warning:** Do not apply the cuff over a wound; otherwise it can cause further injury.
- Do not inflate the cuff on the same limb which other monitoring ME equipment is applied around simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring ME equipment.
- On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure (cuff pressure > 300mmHg or constant pressure > 15mmHg for more than 3 minutes) applied to the arm may lead to an ecchymosis.
- Please check that operation of the device does not result in prolonged impairment of patient blood circulation.
- When measurement, please avoid compression or restriction of the connection tubing.
- The device cannot be used with HF surgical equipment at the same time.
- This device is contraindicated for any female who may be suspected of, or is pregnant. Besides providing inaccurate readings, the effects of this device on the fetus are unknown.
- Too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.
- This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's arm and fingers will become anaesthetic, swollen and even purple due to a lack of blood.
- Keep the monitor at least 20 centimeters away from the human body (especially the head) when the data transmission is proceeding after measurement.
- To enable the data transmission function, this product should be paired to Bluetooth end at 2.4 GHz.
- When not in use, store the device in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.
- **This device may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.**
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in this booklet.
- The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- The patient is an intended operator.
- The patient can measure data and change batteries under normal circumstances and maintain the device and its accessories according to the user manual.
- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- The blood pressure monitor, adaptor and the cuff are suitable for use within the patient environment. If you are allergic to polyester, nylon or plastic, please don't use this device.
- If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the START/STOP button to release the air immediately from the cuff. Loosen the cuff and remove it from your arm.
- If the cuff pressure reaches 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressures reaches 40 kPa (300 mmHg), detach the cuff from the arm and press the START/STOP button to stop inflation.

- Before use, make sure the device functions safely and is in proper working condition. Check the device, do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Do not wash the cuff in a washing machine or dishwasher!
- The service life of the cuff may vary by the frequency of washing, skin condition, and storage state. The typical service life is 10000 times.
- The operator shall not touch output of batteries and the patient simultaneously.
- Don't open or repair the device by yourself in the event of malfunctions. The device must only be serviced, repaired and opened by individuals at authorized sales/service centers.
- Keep the unit out of reach of infants, young children or pets to avoid inhalation or swallowing of small parts. It is dangerous or even fatal.
- Be careful to strangulation due to cables and hoses, particularly due to excessive length.
- At least 30 min required for ME equipment to warm from the minimum storage temperature between uses until it is ready for intended use. At least 30 min required for ME equipment to cool from the maximum storage temperature between uses until it is ready for intended use.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance d away from the equipment.
- There is no luer lock connectors are used in the construction of tubing, there is a possibility that they might be inadvertently connected to intravascular fluid systems, allowing air to be pumped into a blood vessel.
- This product has been factory calibrated. Generally we recommend to calibrate appliance every 2 years if used regularly.

Wireless connectivity and compatibility

- HP TRONIC Zlín, spol. s r.o. can in no way be held liable for any accidental, indirect or subsequent damage, data loss or any other loss incurred as a result of information leakage during communication through BLUETOOTH, Wifi technology.
- 100% successful wireless communication cannot be guaranteed with all Smart devices. Devices have to be compatible with BLUETOOTH specifications set forth by the company Bluetooth SIG, Inc., and their compliance with the specifications has to be certified. However, given the variety of devices available in the market, even if a device is in compliance with the specifications, there might be cases when the character or specifications of the given Smart device do not allow connectivity, or hinder connectivity in any way. This is not a reason to make a warranty claim.
- Radio waves transmitted by devices may interfere with operation of medical equipment. Such interference may cause a failure. Therefore, do not use the device in the following places:
 - in hospitals, trains, planes, petrol stations and in areas where flammable gases may be present.
 - in the vicinity of automatic doors or fire detectors.

The App and the web platform are owned and operated by a third entity (hereinafter referred to as the “Entity”) other than company HP TRONIC Zlín, s r.o. (hereinafter referred to as the “Seller”). Personal data are not transferred between the Entity and the Seller. When signing in to the App, the customer submits their personal data directly to the Entity, which becomes the personal data administrator. The App and related personal data protection are governed by the rules of the Entity.

The company HP TRONIC Zlín, spol. s r.o. uses the services of the company MedM Inc., offering its professional multiplatform apps for the manufacturers of different medical devices and measuring instruments for domestic use. Besides the apps, MedM Inc. offers also a web platform (<https://health.medm.com/>), allowing the customer after logging in (under the same user name and password as in the App) to manage their own measuring in an easier and more comprehensive way (including export, sharing with family members or physician, setting individual limits for measurement results etc.). The App is regularly updated by the company MedM Inc., however, HP TRONIC Zlín, spol. s r.o. has only limited possibilities of its adjustment. The development and all changes in the App are made exclusively by the owner of the App, i.e. by MedM Inc. The company HP TRONIC Zlín, spol. s r.o. has no access to any data related to using the App.

Web platform, available on <https://health.medm.com/>, is only a supplementary service and an option of broadening the use of the product. It is not a part of or an accessory to the ETA4297 pressure gauge! Further, the company HP TRONIC Zlín, spol. s r.o. does not guarantee the functionality of this platform or its localisation into the official language of the customer.

II. GENERAL INFORMATION

You have bought a device, at a first glance looking and acting as an ordinary pressure gauge. However, if connected via bluetooth to any Smart device with Android or iOS operating system with MedM Blood Pressure App installed, the pressure gauge becomes a multifunctional modern assistant to help monitor and analyse the pressure and heart rate.

After linking the pressure gauge with a Smart device, the measurement results are transferred via bluetooth to the **MedM Blood Pressure App**.

App features

- Automatic transfer of measured values to the App;
- Possibility to input the measured values manually;
- Possibility to link the details with individual measurement (e.g. feeling, body position, notes etc.);
- Comprehensive graphical representation of the results;
- Possibility to export the results for selected period to different formats (general CSV, Excel CSV and for the web platform also to PDF, HTML or Excel XLSX)
- The data are automatically saved and backed up in a cloud (on the backed-up servers of the App operator) with an access via own account;
- Possibility to use the App quickly for the visitors even without signing in (so-called visitor mode);
- Memory for a big amount of the saved profiles;
- Allows to switch between blood pressure measurement (European ESH/ESC or U.S. JNC // ACC/AHA);
- Allows the data synchronisation with the GOOGLE FIT (Android) or Health (iOS) apps;
- Designed for Android and iOS (can be downloaded for free from Apple Store / Google Play); Compatible version: Android 4.4 >; iOS 10 >.
- The data used to log in to the App account can also be used for the access to the web platform, adding other options of checking and monitoring your biometrical data.

Note


The pressure gauge can also be fully used without the App as a common pressure gauge in its category. However, using the App enhances the user comfort and the level of checking the blood pressure.

Measurement Principle

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a “zero pressure” equivalent to the atmospheric pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.


What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.

What is the standard blood pressure classification?

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

Blood Pressure (mm Hg)	Optimal	Normal	High-normal	Mild	Moderate	Severe
SYS	<120	120-129	130-139	140-159	160-179	≥180
DIA	<80	80-84	85-89	90-99	100-109	≥110

CAUTION

Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

Irregular Heartbeat Detector

An irregular heartbeat is detected when a heartbeat rhythm varies while the device is measuring systolic pressure and diastolic pressure. During each measurement, blood pressure monitor will keep a record of all the pulse intervals and calculate the average value of them. If there are two or more pulse intervals, the difference between each interval and the average is more than the average value of $\pm 25\%$, or there are four or more pulse intervals, the difference between each interval and the average is more than the average value of $\pm 15\%$, then the irregular heartbeat symbol will appear on the display with the measurement result.

CAUTION

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heart-beat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

III. PREPARATION FOR USE

Remove all the packing material and remove the pressure gauge. Remove adhesive film, stickers, or paper. Choose your power supply - you can use battery power or AC power using an adapter.

Replacing batteries

Open the battery cover on the bottom of the device. Insert a battery with the right polarity and close the cover. Use AAA batteries (4 pcs) for this appliance. If you are not using the appliance for some time, take the batteries out!


Note

- Do not use different type of batteries: Do not use used and new batteries together.
- Do not dispose the batteries in fire. Batteries may explode or leak.
- Remove batteries if the device is not likely to be used for some time.
- Worn batteries are harmful to the environment. Do not dispose with daily garbage.
- Remove the old batteries from the device following your local recycling guidelines.

Descriptions (obr.1)**A – monitor**

- A1 – START/STOP button
- A2 – SET button
- A3 – MEM button
- A4 – LCD display

- A5 – AIR CONNECTOR PLUG
- A6 – BATTERY COMPARTMENT
- A7 – Hole for adapter connection
- A8 – Key lock

B – cuff**C – air hose****DISPLAY DESCRIPTION**

SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic blood pressure	High blood pressure
DIA	Diastolic blood pressure	Low blood pressure
Pul/min	Pulse display	Pulse in beats per minute
kPa	kPa	Measurement Unit of the blood pressure
mmHg	mmHg	Measurement Unit of the blood pressure
	Low battery	Batteries are low and need to be replaced
	Irregular heartbeat	Blood pressure monitor is detecting an irregular heartbeat during measurement.
	Blood pressure level indicator	Indicate the blood pressure level
	Current Time	Month/Day/Year, Hour/Minute
	Heartbeat	Blood pressure monitor is detecting a heartbeat during measurement.
	User 1/User 2	Start measurement for User 1/User 2
	The average value	The average value of the latest three records
	Bluetooth icon	The bluetooth icon blinks when the Bluetooth is working
	Motion indicator	Motion may result in an inaccurate measurement
	Memory	Indicate it is in the memory mode and which group of memory it is.
	Bluetooth connection	It blinks when the bluetooth connection fails or the data is not sent.

IV. INSTRUCTIONS FOR USE

A) SETTING DATE, TIME AND MEASUREMENT UNIT

It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory. (The setting range of the year: 2018-2058 time format: 24H)

- 1) When the monitor is off, hold pressing “**SET**” button to enter the mode for year setting.
- 2) Press the “**SET**” to change the [YEAR]. Each press will increase the numeral by one in a cycling manner.
- 3) When you get the right year, press “**SET**” to set down and turn to next step. Repeat step 2 and 3 to set the [MONTH] and [DAY].
- 4) Repeat step 2 and 3 to set the [HOUR] and [MINUTE].
- 5) Repeat step 2 and 3 to set the [UNIT].
- 6) After the unit is setted, the LCD will display “**dONE**” first, then display all the settings you have done and then it will turn off.

B) SELECT THE USER

- 1) When the monitor is off, long press the „MEM“ button to enter user setting mode.
- 2) Then press „MEM“ button again, select the user ID between user 1 and user 2.
- 3) After selecting the suitable user ID, press „SET“ button to confirm.
Then the LCD will turn off.

C) TIE THE CUFF (fig. 2, 3 a 4)

- 1) Remove all jewelry, such as watches and bracelets from your left arm. **Note:** If your doctor has diagnosed you with poor circulation in your left arm, use your right arm.
- 2) Roll or push up your sleeve to expose the skin. Make sure your sleeve is not too tight.
- 3) Hold your arm with your palm facing up and tie the cuff on your upper arm, then position the tube off-center toward the inner side of arm in line with the little finger. Or position the artery mark Φ over the main artery (on the inside of your arm). Fig.2
Note: Locate the main artery by pressing with 2 fingers approximately 2 cm above the bend of your elbow on the inside of your left arm. Identify where the pulse can be felt the strongest. This is your main artery. Fig.3
- 4) The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.
- 5) Sit comfortably with your tested arm resting on a flat surface. Place your elbow on a table so that the cuff is at the same level as your heart. Turn your palm upwards. Sit upright in a chair, and take 5-6 deep breaths. Fig.4
- 6) Helpful tips for Patients, especially for Patients with Hypertension:
 - Rest for 5 minutes before first measuring.
 - Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
 - Take the measurement in a silent room.
 - The patient must relax as much as possible and do not move and talk during the measurement procedure.
 - The cuff should maintain at the same level as the right atrium of the heart.
 - Please sit comfortably. Do not cross your legs and keep your feet flat on the ground.
 - Keep your back against the backrest of the chair.

For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

D) START THE MEASUREMENT

- 1) When the monitor is off, press the “START/STOP” to turn on the monitor.


Adjust the zero

Inflating and measuring

Display and save the results

2) Press the “**START/STOP**” to power off, otherwise it will turn off within 1 minute.

Note

To avoid unwanted starting, long press the **A8** button to activate the button lock.
Long press again to cancel the button lock.

E) LINKING THE PRESSURE GAUGE WITH A SMART DEVICE

First, download the application named **MedM Blood Pressure** to your smart device. It is available on Google Play for Android, and in the App Store for iOS.


Basic summary

First, after downloading the app, select the mode of use (with or without signing in; <1>). You can use the app without signing in, however, the data are saved only to your smart device and if the app is de-installed, the data are lost. If signed in via e-mail, the data are saved on the server, so you have the access to your data from any other devices, including the access to the web platform.

For signing in you must enter your e-mail, password and agree with the Terms and Conditions of the app provider.

<2>. After activation of the app, you can set your PIN code, without which the app cannot be activated. Before you start to use the app, select in the Settings the <3> gauge, or the methodology of blood pressure measurement <4, 5>. Allow the app to get the accesses required (e.g. storage site, bluetooth etc.)

Menu

Current measurement

Last measurement date and time

Brief overview of last measurements

Detailed overview of last measurements

Selection and administration of users

Manual input of the measured values

Details and comments to the measurement

Blood pressure values graph

Heart rate values graph

30 days report view

The main display shows a blood pressure reading of 109 / 68 64. The date and time of the last measurement is 4:44 13 Feb 2020. Below this is a list of recent measurements:

Date	Time	SYS	DIA	HR
13 Feb	14:44	109	68	64
13 Feb	14:34	108	69	66
13 Feb	14:28	131	71	72
13 Feb	14:27	152	76	69

The 30-day report shows a 54% optimal blood pressure rate for the period 14 Jan - 13 Feb. The report includes a bar chart showing the number of measurements in different categories: Optimal (54%), Normal (15%), and High Normal (23%).

Linking procedure

- 1) Activate the Bluetooth and then the **MedM** app.
- 2) In the section Menu – My devices click on “**Add new**”.
- 3) The app will start searching for an active device in its surroundings. In this moment activate on the switched off pressure gauge the linking function by long pressing the button **A3 START/STOP**.
- 4) The found device is displayed, click on it and confirm the selection “add to the device” in the lower part of the display.
- 5) If linking is successful, the pressure gauge is displayed “**⌘ + []**” including user ID. If linking fails, the pressure gauge is displayed “**⌘ + []**” including user ID.
- 6) When the linking process is completed, the pressure gauge is automatically switched off.

TIP

- Before confirming the addition of the device, choose profile for which the linking will be active.
- The app is intuitive and allows to monitor the measurement results in different modes. It will help you better monitor your pressure and heart rate. You can show the results in the app (or exported and printed) e.g. to your physician.

Automatic data transfer


In the section **Menu – My devices - (linked devices)** - you can click on the setting icon in the right upper part and de/activate the function of an automatic transfer of the measured values. If the function is active, the measured data is immediately transferred to the active app.

If the function is not active, the results must be transferred manually by clicking on , then on the pressure gauge icon and perform other measurement. When the measurement is completed, the pressure gauge will be in the mode of the result transfer and the transfer of selected results can be selected and confirmed.

In the same part of the setting you can select the user profile, for which the linking will be valid (profile 1, profile 2, or both profiles).

Data export

In the section **Menu - Export** there are options of the data export to the CSV files. For more export options (PDF, HTML či Excel XLSX) use the access via the web platform.

Internet platform (<https://health.medm.com>)

The company MedM provides the manufacturers with its apps designed for different product types (pressure gauges, thermometers, weight scales, glucometers etc.). If a customer owns more supported devices and uses different types of MedM applications for them, they can also use our internet platform, in which they can administer the data from all these devices in one place. Although you own only a pressure gauge, you can also use the benefits of the platform.

In principle, you can find here the same data as shown in the app + some extra benefits:



- Better orientation in the results and graphs;
- Option of setting/adjustment of individual graph limits (i.e. what results are suitable for you and what represent high/low blood pressure);
- More data export options;
- Option of setting a notification (e.g. setting up an e-mail with a new measurement, or if the measurement exceeds the set limits);
- Option to share the data with other persons (family members, physician etc.); the function can be used e.g. to have the results (or only exceeding the defined values) sent to the grandchildren/ children of the elderly, so they have a better awareness of the health condition of their family members (note: a person who wants to have an access to the data of other users, must have their own account created)

For logging in, the same account can be used as the account used for logging in to the app **MedM Blood Pressure**.


F) RECALL THE RECORDS

- 1) When the monitor is off, please press the “MEM” to show the average value of the latest three records. If the records are less than three groups, it will display the latest record first.
- 2) Press the “MEM” or “SET” to get the record you want.

The date and time of the record will be shown alternately.

CAUTION

The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (60) is dropped.

– To view the data of the other user, simply switch the user profile.

G) DELETE THE RECORDS

If you did not get the correct measurement, you can delete all results by following steps below.

- 1) Hold pressing “MEM” for 3 seconds when the monitor is in the memory recall mode, the flash display will show „dEL ALL“.
- 2) Press “SET” to confirm deleting, the monitor will display “dEL done” and then turn off.
- 3) If you don’t want to delete the records, press “START/STOP” to escape.
- 4) If there is no record, the display will show „ - - - “.

V. TROUBLESHOOTING OF POSSIBLE PROBLEMS

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the products not operating as you think it should, check here before arranging for servicing.

PROBLEM	SYMPTOM	CHECK THIS	REMEDY
Error message	E 01 shows	The cuff is too tight or too loose.	Refasten the cuff and then measure again.
	E 02 shows	The monitor detected motion while measuring.	Movement can affect the measurement. Relax for a moment and then measure again.
	E 03 shows	The measurement process does not detect the pulse signal.	Loosen the clothing on the wrist and then measure again.
	E 04 shows	The treatment of the measurement failed.	Relax for a moment and then measure again.
	EExx, shows on the display.	A calibration error occurred.	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty for contact information and return.
Warning message	“out ” shows	Out of measurement range	Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your physician.

VI. MAINTENANCE

Dust environment may affect the performance of the unit. Please use the soft cloth to clean the whole unit before and after use. Don't use any abrasive or volatile cleaners.

In order to get the best performance, please follow the instructions below.

- Put in a dry (including adaptor) place and avoid the sunshine
- Avoid touching water, clean it with a dry cloth in case.
- Avoid intense shaking and collisions
- Avoid dusty and unstable temperature environment
- Using wet cloths to remove dirt
- Do not attempt to clean the reusable cuff with water and never immerse the cuff in water.

More extensive maintenance or maintenance requiring intervention in the inner parts of the appliance must be carried out by a professional service! Failure to comply with the instructions of the manufacturer will lead to expiration of the right to guarantee repair! Hereby, ETA a.s. declares that the radio equipment type ETA4297 is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at eta webpage.

VII. ENVIRONMENTAL PROTECTION

If the dimensions allow, there are labels of materials used for production of packing, components and accessories as well as their recycling on all parts.

The symbols specified on the product or in the accompanying documentation mean that the used electrical or electronic products must not be disposed of together with municipal waste. For proper disposal, hand them over at special collection places where they will be accepted free of charge. Suitable disposal of the product can help to maintain valuable natural resources and to prevent possible negative impacts on the environment and human health, which could be the possible consequences of improper waste disposal.

Ask for more details at the local authorities or at a collection site. Fines can be imposed for improper disposal of this type of waste in agreement with the national regulations.

Remove a discharged battery and dispose of it in a suitable and safe way through the special collection networks. Never dispose of the battery by burning!

VIII. TECHNICAL DATA

Parameter	Description
Measurement mode	Oscillographic testing mode; Protection against
Measurement range	Rated cuff pressure: 0 mm Hg ~ 299 mmHg (0 kPa ~ 39.9 kPa) Measurement pressure: SYS: 60 mm Hg ~ 230 mmHg (8.0 kPa ~ 30.7 kPa) DIA: 40 mmHg ~ 130 mmHg (5.3 kPa ~ 17.3 kPa) Pulse value: (40-199) beat/minute
Accuracy	Pressure: 5 °C - 40 °C within ±0.4 kPa (3 mmHg) Pulse value: ±5 %
Normal working condition	A temperature range of: +5 °C to +40 °C A relative humidity range of: 15 % to 90 % Non-condensing, but not requiring a water vapour partial pressure greater than 50 hPa An atmospheric pressure range of: 700 hPa to 1060 hPa
Storage temperature	-20 °C ~ +60 °C, less than 93 % HR non-condensing, at a water vapour pressure up to 50 hPa
Measurement perimeter of the upper arm	About 22 cm - 42 cm
Protection against ingress of water	IP21 It means the device could protected against solid foreign objects of 12.5mm and greater, and protect against vertically falling water drops
Power supply	4 x 1,5 V AAA (included); Adaptor (model BLJ06W060100P1-U); AC 100-240V; 50/60Hz 0.2 A Max; Napětí 6 V \approx 1000 mA
Net weight of device	0,25 kg
Size of the product (mm)	118 x 107 x 103
Degree of protection	Type BF applied part

Parameter	Description
Software version	A01
Bluetooth Module No.:	LS51802
RF Frequency Range:	2402 MHz to 2480 MHz
Output Power Range:	≤4dBm
Supply Voltage:	2 V - 3.6 V
Transmitting Distance:	10 meters

WARNINGS AND SYMBOLS USED ON THE APPLIANCE, PACKAGING OR IN THE INSTRUCTIONS MANUAL:

HOUSEHOLD USE ONLY. DO NOT IMMERSE IN WATER OR OTHER LIQUIDS.

TO AVOID DANGER OF SUFFOCATION, KEEP THIS PLASTIC BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE THIS BAG IN CRIBS, BEDS, CARRIAGES OR PLAYPENS. THIS BAG IS NOT A TOY.

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "TYPE BF APPLIED PARTS"
CE 0123	Symbol for "COMPLIES WITH MDD 93/42/EEC REQUIREMENTS"		Symbol for "ENVIRONMENT PROTECTION - Electrical waste products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "MANUFACTURER"		
	Symbol for "SERIAL NUMBER"		
	Symbol for "DIRECT CURRENT"		Symbol for "Authorised Representative in the European Community"
	Symbol for "MANUFACTURE DATE"		Symbol for "RECYCLE"
	Caution: These notes must be observed to prevent any damage to the device.		The Green Dot is the license symbol of a European network of industry-funded systems for recycling the packaging materials of consumer goods.
F1	T1A/250V Φ3.6*10CCC		Symbol for "Class II Equipment"
	For indoor use only		

Device Classification

Battery Powered Mode: Internally Powered ME Equipment

AC Adaptor Powered Mode: Class II ME Equipment

Manufactured by: Guangdong Transtek Medical Electronics Co., Ltd.
Zone A, No.105 ,Dongli Road, Torch Development District
Zhongshan,528437,Guangdong,China

Authorized European Representative:

MDSS - Medical Device Safety Service GmbH; Schiffgraben 41, 30175 Hannover,
Germany

Importer: ETA a.s.; Křižíkova 148/34; Praha 8 - Karlín, 186 00

Complied Standards List

Risk management	EN ISO 14971:2012 / ISO 14971:2007 Medical devices - Application of risk management to medical devices
Labeling	EN ISO 15223-1:2016 / ISO 15223-1:2016 Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 1 : General requirements
User manual	EN 1041:2008 +A1:2013 Information supplied by the manufacturer of medical devices
General Requirements for Safety	EN 60601-1:2006+A1:2013/ IEC 60601-1:2005+A1:2012 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance EN 60601-1-11:2015/ IEC 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Electromagnetic compatibility	EN 60601-1-2:2015/ IEC 60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests
Performance requirements	EN ISO 81060-1:2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type EN 1060-3:1997+A2:2009 Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems IEC 80601-2-30:2009+A1:2013 Medical electrical equipment- Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers
Clinical investigation	EN 1060-4:2004 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers ISO 81060-2:2013 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated measurement type
Usability	EN 60601-1-6:2010+A1:2015/IEC 60601-1-6:2010+A1:2013 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability IEC 62366-1:2015 Medical devices - Part 1: Application of usability engineering to medical devices
Software life-cycle processes	EN 62304:2006/AC: 2008 / IEC 62304: 2006+A1:2015 Medical device software - Software life-cycle processes
Bio-compatibility	ISO 10993-1:2009 Biological evaluation of medical devices- Part 1: Evaluation and testing within a risk management process ISO 10993-5:2009 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity ISO 10993-10:2010 Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization

EMC GUIDANCE

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment TMB-1583-BS, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

- 1) All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
- 2) Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/flicker emissions IEC 61000-3-3	Comply

Table 2

Guidance and manufacturer's declaration – electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines 1 kV signal input/output 100 kHz repetition frequency	±2 kV for power supply lines 1 kV signal input/output 100 kHz repetition frequency
Surge IEC61000-4-5	±0.5 kV, ±1 kV differential mode 0.5 kV, ±1 kV,±2 kV common mode	±0.5 kV, ±1 kV differential mode 0.5 kV, ±1 kV,±2 kV common mode
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°.0 % UT; 250/300 cycle	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°.0 % UT; 250/300 cycle

Guidance and manufacturer's declaration – electromagnetic Immunity

Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V; 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V; 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz; 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz; 80 % AM at 1 kHz

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380-390	TETRA 400	Pulse modulation b) 18Hz	1,8	0,3
450		430-470	GMRS 460, FRS 460	FM c) \pm 5kHz deviation 1kHz sine	2	0,3	28
710		704-787	LTE Band 13,17	Pulse modulation b) 217Hz	0,2	0,3	9
745							
780							
810		800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18Hz	2	0,3	28
870							
930							
1720		1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4,25; UMTS	Pulse modulation b) 217Hz	2	0,3	28
1845							
1970							
2450		2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
5240		5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9
5500							
5785							

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© DATE 16/3/2020